



# 25mL Centrifuge Tubes

Closing the obvious gap between 15mL and 50mL tubes, they are ideal for any application that requires more than 15mL and less than 50mL, including midpreps, general sample prep, centrifugation, mixing and storage.

## **GenFollower Biotech CO., LTD.**

Address: No.273, Yuexiuzhong road, cao'e street, Shangyu area, Shaoxing, Zhejiang, China  
Phone: +86-575-82936953 / Fax: +86-575-82032762  
Email: [GenFollower@GenFollower.com](mailto:GenFollower@GenFollower.com)



- ⊕ **Maximum force up to 12,000xg**  
> High clarity polypropylene can withstand up to 12,000 x g maximum RCF

- ⊕ **Lower Height and large diameter design**

> Better Sample Availability & Less Contamination : the insertion depth of the pipette into the tube is much shorter compared to 15 mL or 50 mL conical tubes, the risk of cross contamination between the tube and the pipetting system are thus reduced to a minimum.

> Low-profile design that minimizes shipping, storage and waste.

> Reduced height and large diameter design (which is identical to 50mL tubes) allows for faster re-suspension of pellets and 100% recovery of small samples.

- ⊕ **Working Temperature: -86°C to 100°C**

- ⊕ **Sterile: RNase, DNase and Pyrogen Free**

> Printed graduations from 5ml to 25ml & with large writing area



| Size | Packaging       | Color | Class  | Cat.#   | Quantity per Case                      | Size/Case(CM) |    |    | Weight/Case(Kg) |
|------|-----------------|-------|--------|---------|--|---------------|----|----|-----------------|
|      |                 |       |        |         |  | L             | w  | H  |                 |
| 25mL | 25PCS/Bag       | Blue  | Common | CTB25BL | 8 Bags<br>(25 PCS/Bag, 8 Bags/Case)    | 43            | 35 | 21 | 2.5             |
| 25mL | 25PCS/Foam Rack | Blue  | Common | CTF25BL | 8 Racks<br>(25 PCS/Rack, 8 Racks/Case) | 43            | 35 | 21 | 2.7             |